

## Competent cells (DH5 $\alpha$ and BL21).v1

This protocol is modified based on NAR 16(8) 1988 page 3580

### Procedures

#### **Preparation**

1. Grow an overnight culture of DH5 $\alpha$  or BL21 bacteria on an LB plate (WITHOUT antibiotics).
2. Inoculate the bacteria into 3-ml LB medium in a 13-ml test tube and shake at 250 rpm at 37°C O/N.
3. Pipette 2 ml of the overnight culture into 500 ml of LB in a 2-liter flask. Shake at 220 rpm at 37°C to allow to grow to 0.3-0.6 at OD600. It takes about 3 hours. Need to check periodically, especially when the OD reading is getting close to 0.3-0.6 (the exponential growth period).
4. Spin the culture at 1,000 x g for 10 min at 4°C.
5. Resuspend the pellet gently in 50 ml of ice-cold TSB.
6. Incubate on ice for 10 min
7. Aliquot in 1.5-ml Eppendorf tubes, quick freeze in an ethanol-dry ice bath and store at -80°C.

#### Transformation

1. Thaw the bacteria on ice. Transfer 100  $\mu$ l to each 13-ml tube on ice.
2. Mix with 100 pg of plasmid DNA or 2  $\mu$ l of ligation mixture.
3. Incubate on ice for 20 min. Heat shock for 45 seconds at 43°C. Put on ice for 2 min.
4. Add 0.9 ml of TSBG and shake at 37°C for 1 h.
5. Plate the cells on plate with appropriate antibiotics.

#### Reagents and Solutions

##### **TSB (500 ml)**

LB broth (pH 6.1)	465 ml
PEG (MW 3500)	50 g (10%)
1 M MgCl <sub>2</sub>	5 ml (10 mM)
1 M MgSO <sub>4</sub>	5 ml (10 mM)

Autoclave and add 25 ml DMSO to a final concentration of 5% and store at 4°C.

##### **TSBG (50 ml)**

TSB	50 ml
40% (2.22 M) Glucose	450 $\mu$ l

Store at 4°C.