

**Preparation of Competent Cells (JM109) ---Acid Salt Method****Procedures:**

1. **Streak the desired bacterial strain** on M-9 plates and incubate at 37°C overnight.
2. **Inoculate 25 ml of LB medium (NO antibiotics!)** with a single colony and incubate at 30°C overnight with vigorous shaking.
3. **Inoculate 500 ml of LB medium** in a 3-liter flask with 5ml of cells from the overnight culture and shake the culture at 150-200 rpm at 30°C until the A<sub>600</sub> reaches 0.45-0.55.
4. **Chill the cells** in ice water for 2 hrs and collect by centrifugation in 150 ml Corex tube at 4000 rpm (2,500 x g) with rotor GSA (r=5.75 in) for 15-20 minutes at 4°C.
5. **Resuspend the cells** in 10-20 ml of ice-cold Trituration buffer and then dilute to 500 ml with the same solution and incubate the cells on ice for 45 min.
6. **Centrifuge the cells** at 3,300 rpm (1,800 x g) for 10 min at 4°C and gently resuspend in 50 ml of ice cold Trituration buffer.
7. **Pool the cell and add 80% glycerol dropwise** with gentle swirling to a final concentration of 15% (v/v). Make 0.2-1ml aliquots of the cells, freeze on dry ice and store at -70°C.

**Solutions:****Trituration buffer** (600 ml):

6.66 g	CaCl <sub>2</sub> (MW 110.99)	100 mM (or 8.82 g for CaCl <sub>2</sub> ·2H <sub>2</sub> O)
8.54 g	MgCl <sub>2</sub> ·6H <sub>2</sub> O (FW 203.3)	70mM
3.27 g	NaAc·3H <sub>2</sub> O (FW 136.08)	40mM pH 5.5

This solution should be prepared fresh and filter-sterilized for each use.