# **Direct Competitive ELISA (cytochrome c)**

3:32 PM 9/15/09

REAGENTS

Z. Yan

| <u>10X PBS (1 liter)</u> |      | <u>PBSN</u>      |        |
|--------------------------|------|------------------|--------|
| NaCl                     | 80 g | 1X PBS           | 100 ml |
| KC1                      | 2 g  | NaN <sub>3</sub> | 50 mg  |

1 M H<sub>2</sub>SO<sub>4</sub>

Na<sub>2</sub>HPO<sub>4</sub>·7H<sub>2</sub>O 11.5 gMix well and store at room temperature.

KH<sub>2</sub>PO<sub>4</sub>

Add  $H_2O$  to 1 liter (Working solution pH ~7.3).

| BBS (2 liters)   |      | Blocking Reagent |             |
|--|------|------------------|-------------|
| $H_3BO_4$  | 21 g | BBS              | 500 ml      |
| NaCl   | 14g  | Tween 20         | $250 \mu 1$ |
| Dissolve in H <sub>2</sub> O, pH to 8.5 and volume to 2 liters |      | 0.5 M EDTA       | 1 ml        |
| and store at room temperature.                                 |      | BSA              | 1.25 g      |
|  |      | $NaN_3$          | 0.25 g      |
|  |      | Store at 4°C.    |             |

### **Homogenization Buffer**

1 M KCl 17.5 ml 18 M H<sub>2</sub>SO<sub>4</sub> 22.22 ml Add to H<sub>2</sub>O to final volume of 400 ml. Ingredients for 100 ml of H.Bx for luciferase

and CAT.

### **PROCEDURE**

## **Criss-Cross Serial Dilution Analysis**

- Prepare antigen (Ag) (rat heart cytochrome c, Sigma C-7892) in PBSN by 1:2 serial dilution. Add 50 µl of Ag to each well of a row of a 96-well EIA microplate (Bio-Rad 224-0096). The highest concentrion (10  $\mu$ g/ml) to row A, the second highest concentration to row B and the lowest concentration (0  $\mu$ g/ml) to row H. Incubate at room temperature overnight.
- Rinse the coated plate with double distilled water three times.
- Add 150  $\mu$ l of blocking reagent to each well and incubate at room temperature for 30 min.
- Rinse the plate with double distilled water three times.
- Add 50 µl of blocking reagent to each wells except the wells in the first column.
- Add 100 µl of primary antibody (Ab) (mouse monoclonal Ab against rat cytochrome c, 6H2.B4, IgG 2a, K, 0.7 mg/ml from Dr. R. Jemmerson, 1:100 in blocking reagent) to the wells in the first column and incubate at room temperature for 2 hr.
- Rinse the plate with double distilled water three times.
- Add 50 µl of secondary Ab (goat anti-mouse IgG conjugated to horse radish peroxidase in blocking reagent, Axell #JGM035003 diluted in PBS-Tween, 1:5000 in blocking reagent) and incubate at room temperature for 2 hr.
- 9. Rinse the plate with double distilled water three times.
- 10. Add 75  $\mu$ l of blocking reagent and vortex and incubate at room temperature for 10 min.
- 11. Rinse the plate with double distilled water three times.
- 12. Add 100 µl of TMB substrate (Bio-Rad, 172-1066, mix A solution with B solution at 9:1 before use) and incubate at room temperature for 10 min.
- 13. Add 100 µl of 1 M H<sub>2</sub>SO<sub>4</sub> and read absorbance at 450 nm in a microplate reader (Bio-Rad) using the mean of the row H values as zero.
- 14. Plot the data for each Ag dilution (A450 for Y-axis and log1/D for X-axis). The concentrations of Ag that gives a maximal A450 value would be appropriate for coating the plate. The dilution of primary Ab that

# **Direct Competitive ELISA (cytochrome c)**

3:32 PM 9/15/09

Z. Yan

gives the 1/2 of the maximal reading at A450 would be used for sensitive detection of cytochrome c in the direct competitive ELISA.

## Microplate coating

- 1. Add 400  $\mu$ l of cytochrome c (500  $\mu$ g/ml in H.Bx. for luciferase and CAT) to 50 ml of PBSN and the final concentration for cytochrome c is 4  $\mu$ g/ml.
- 2. Add 50  $\mu$ l of the Ag in PBSN to each well of a microplate and incubate at room temperature overnight. The coated plate can be stored at 4°C for several months.

## **Muscle homogenization**

- 1. Weigh out ~100 mg frozen muscle while frozen and put into a 15 ml polyethylene tube (Sarstedt 60.541).
- 2. Add 10 ml/g 1X homogenization buffer (H.Bx for luciferase and CAT with 175 mM KCl).
- 3. When the muscle is thawed on ice, homogenize the sample at setting 7 with a Polytron homogenizer for 10" x 3 with 10" intervals in ice-water bath.
- 4. Freeze-thaw the samples three times using -80°C freezer or liquid nitrogen and ice.
- 5. Centrifuge the samples at 11,000 rpm (15,000 x g) for 15' at 4°C.
- 6. Transfer the supernatant to an eppendorf tube and store at -20°C.

### **ELISA**

- 1. Rinse the coated plate with double distilled water three times.
- 2. Add 150  $\mu$ l of blocking reagent to each well and incubate at room temperature for 30 min.
- 3. Rinse the plate with double distilled water three times.
- 4. Prepare cytochrome c standard in homogenization buffer (100 ng/ $\mu$ l is the highest concentration and 2:3 series dilution for the rest 9 tubes and homogenization buffer as the 0 ng/ $\mu$ l control).
- 5. To wells in a NEW microplate, add 100  $\mu$ l of blocking reagent containing 1:9000 primary Ab and add 11.11  $\mu$ l of standard or sample in duplicates to each well. Incubate at room temperature for 30 min.
- 6. To coated well add 50  $\mu$ l of the above mixture or 50  $\mu$ l of blocking reagent (as negative control to zero the reading) and incubate at room temperature for 2 hr.
- 7. Rinse the plate with double distilled water three times.
- 8. Add 50  $\mu$ l of secondary Ab and incubate at room temperature for 2 hr.
- 9. Rinse the plate with double distilled water three times.
- 10. Add 75 µl of blocking reagent and vortex and incubate at room temperature for 10 min.
- 11. Rinse the plate with double distilled water three times.
- 12. Add 100  $\mu$ l of TMB substrate and incubate at room temperature for 10 min.
- 13. Add 100 µl of 1 M H<sub>2</sub>SO<sub>4</sub> and read absorbance at 450 nm in a microplate reader.

### REFERENCES

- 1 Frederick M. Ausubel, Roger Brent, Robert E. Kingston, David D. Moore, J. G. Seidman, John A. Smith, and Kevin Struhl. *Current Protocols in Molecular Biology*, Chapter 11: Immunology. John Wiley & Sons, Inc., 1994.
- 2. Jemmerson, R., C. Mueller, and D. Flaa. Differences in heavy chain amino acid sequences affecting the specificity of antibodies for variants of cytochrome c. *Molecular Immunology*. 30:1107-1114, 1993.
- 3. Holloszy, J. O. Biochemical adaptations in muscle. Effects of exercise on mitochondrial oxygen uptake and respiratory enzyme activity in skeletal muscle. Journal of Biological Chemistry. 242:2278-2282, 1967.