

Fixation for immunocytochemistry for tissue samples.v3

This fixation procedure will allow you to process your samples to a point that they can be store away at -20°C temperature before you decide whether you want to proceed immunocytochemical labeling.

Procedures:

1. **Contact Jan A. Redick** at the Advanced Microscopy Facility (Tel: (434) 924-2524, email: jar@virginia.edu) to set up the date and time for using the 45°C oven.
2. Euthanize the animals **under anesthesia** with cervical dislocation.
3. Dissect and harvest the tissue **as quickly** as possible.
4. Put a few drops of fixative (see below) on wax or plastic petri dish and immediately **put the harvested tissue into the fixative**. Then cut the tissue specimen into small pieces (1 mm cubes or smaller). For skeletal muscle, it is ideal to have a thin longitudinal strip.

THE FOLLOWING PROCEDURES ARE DONE AT 4°C

5. Immediately **put the sample into 4 ml of ice-cold 4% paraformaldehyde + 0.2% glutaraldehyde in 1X PBS (pH 7.2)**, in a 5-ml polypropylene tube on ice and gently rotate at 4°C for 1 hour.
6. Wash samples in 2 ml of ice-cold **1X PBS** 3 x 10 min with rocking.
7. Wash samples in 2 ml of ice-cold **40% ETOH** for 10 min with rocking.
8. Wash samples in 2 ml of ice-cold **60% ETOH** for 10 min with rocking.
9. Wash samples in 2 ml of ice-cold **80% ETOH** for 10 min with rocking.
10. Wash samples in 2 ml of ice-cold **100% ETOH** for 2 x 10 min with rocking.

THE FOLLOWING PROCEDURES ARE DONE AT ROOM TEMPERATURE

11. Incubate (with rocking) the samples in 2 ml of **2:1 100% ETOH:LRW** for ½ -1 hour.
12. Incubate (with rocking) the samples in 2 ml of **1:1 100% ETOH:LRW** for 2-3 hrs (or overnight) (OK stopping point for convenience)
13. Incubate (with rocking) the samples in 2 ml of **1:2 100% ETOH:LRW** for 2 hrs (or overnight) (OK stopping point for convenience)
14. Incubate (with rocking) the samples in 2 ml of **1:4 100% ETOH:LRW** for 2 hrs (or overnight) (OK stopping point for convenience)
15. Incubate (with rocking) the samples in 2 ml of **100% LRW** for 1-2 hrs (or overnight) (OK stopping point for convenience)
16. Incubate (with rocking) the samples in 2 ml of **100% LRW** for 2-3 hrs.
17. Incubate (with rocking) the samples in 2 ml of **100% LRW** for 2-3 hrs.
18. Embed the samples in **gelatin capsules** (not in BEEM capsules or Eppendorf tubes): First, write labels in pencil (not pen or printer ink), and put labels inside capsule, near the top. Fill capsule to top with **fresh LRW**. Add tissue and let sink to bottom. Push capsule tops on snugly (NOTE: Oxygen will inhibit polymerization).
19. Polymerize the samples in capsule for 24 hours in a **45°C oven**. (Note: do not polymerize at the usual 60°C...too hard on antigens)

Solution preparation:

8% paraformaldehyde in 1X PBS

1. Wear gloves, goat and goggles.

2. Add 70 ml of ddH₂O into a 250-ml beaker.
3. Add 8 g paraformaldehyde and mix them on a stir plate in a fume hood.
4. Add NaOH (~13-15 pellets) to bring pH to ~12 and wait until the paraformaldehyde powder gets completely dissolved.
5. Add 10 ml of 10X PBS, pH 7.2.
6. Add 10N HCl to bring pH back to ~8.0 and then use 1N HCl to bring it to 7.2.
7. Double check pH to make sure it is 7.2 (confirm with pH paper)
8. Add about 10 ml ddH₂O to a final volume of 100 ml.
9. Filter sterilize the solution and store at 4°C up to 1 month.

0.4% Glutaraldehyde in 1X PBS

1. Add 5 ml of 10X PBS, pH 7.2, to 50 ml conical tube.
 2. Add 44.2 ml of ddH₂O.
 3. Add 0.8 ml of 25% glutaraldehyde* to a final volume of 50 ml.
 4. Make this solution fresh for each use.
- * Pippette up and down the glutaraldehyde to ensure the solution is homogenous before adding.

4% paraformaldehyde, 0.2% glutaraldehyde in 1X PBS

Mix 8% paraformaldehyde in 1X PBS with 0.4% glutaraldehyde in 1X PBS 1:1 to make the final solution. This solution should be made fresh for each use. Label it “**For immunogold**”.

ETOH:LWR mix

Make this fresh before use.