

Protocol of serum or plasma BUN assay.v2

This is a modified assay for quantification of blood urine nitrogen according to DetectX® BUN Colorimetric Detection Kit (K024-H1) Product Protocol.

Sample Preparation

1. Prepare mouse serum or plasma according to protocols described on our website.
2. Dilute serum and plasma in ddH₂O \geq 1:10 and \geq 1:20, respectively (recommended).
3. Dilute serum: Pipet 4 μ l serum in 76 μ l of ddH₂O in our experience (1:20).

Standard Preparation

Use distilled or deionized water to prepare standard from BUN standard stock (100 mg/dl)

| | Std 1 | Std 2 | Std 3 | Std 4 |
|-------------------------------|--------------------|------------------|--------------------|-------|
| ddH ₂ O (μ l) | 138.7 | 37.5 | 37.5 | 75 |
| Addition | 11.3 μ l Stock | 75 μ l Std 1 | 37.5 μ l Std 2 | N/A |
| Final conc (mg/dl) | 7.5 | 5 | 2.5 | 0 |

※ Use the standards within 2 hours of preparation.

Procedure

1. Run all samples and standards in duplicates.
2. Pipet 33 μ l of samples or appropriate standards in a 96 well plate in duplicates.
3. Add 50 μ l of Color Reagent A to each well using a repeater pipette.
4. Add 50 μ l of Color Reagent B to each well using a repeater pipette.
5. Incubate at room temperature for 30 minutes.
6. Read the optical density at 450 nm.

Calculation

Average the duplicate OD readings for each standard and sample. Create a standard curve after subtracting the mean OD's for the blank (Std 4). Calculate sample concentrations using the standard curve and multiplying by the dilution factor.