

SDH staining protocol. v1

This protocol is based on J. Histochem. (Nachals, NM et al. 5. 420-430, 1957) for staining skeletal muscle cross-section to determine if there is a change of mitochondrial SDH activity as an index of oxidative phenotype,

Procedure

1. **Sacrifice mice** by cervical dislocation under anesthesia followed by muscle sample harvesting. Cut the muscle across at the belly.
2. Freeze the muscle in a position that the muscle fibers are perpendicular to the bottom of the mold in OCT freezing compound (Tissue-Tek, Sakura) pre-chilled in liquid nitrogen. Wrap the tissue block in an aluminum foil and store in -80°C.
3. Cut skeletal muscle cross-sections (5 μ m) using a cryostat and mount it on a **positive charged glass slide**, wrap in an aluminum foil and store at -80°C until use.
4. Place the slide with muscle section in a chamber, and incubate in prewarmed (37C°) **Incubation Solution** for 40 min at 37C° in a water bath. This incubation time can be extended to 60 min, if the staining is insufficient.
5. Wash the samples in distilled water (1 min x 3 times).
6. Wipe off water from slide glass (**DO NOT touch muscle section**), mount with **Glycerol based mounting medium**, and seal with nail polish.

Reagent

Incubation solution (20 ml)

Reagent	Source/Cat#	MW/Conc.	Quantity	Final Conc.
Sodium succinate	Sigma S2378	270.14/0.2 M	5 ml	50 mM
Phosphate buffer		0.2 M	5 ml	50 mM
Nitroblue tetrazoliumand	Sigma N5514	1 mg/ml	10 ml	0.5 mg/ml

Phosphate buffer stock solution (100 ml)

Reagent	Source/Cat#	MW/Conc.	Quantity	Final Conc.
KH ₂ PO ₄	Fisher P380	136.09/1 M	12ml	0.12M
Na ₂ HPO ₄	Fisher S374	141.96/1 M	88ml	0.88M

Glycerol based mounting medium (50 ml)

Reagent	Source/Cat#	MW/Conc.	Quantity	Final Conc.
Gelatine	Sigma G1890		2.5 g	5%
Glycerol	Fisher, G33-1		25 ml	50%
Distilled H ₂ O			25ml	