Streptozocin Injection.v1

Procedures:

- 1. IMPORTANT: Inform Vivarium supervisor (Laurie Edwards at <a href="https://linear.com/ltm.co
- 2. Remove the food from mice (DO NOT remove the waster bottle) and fast the mice overnight starting around 5 pm (~16 hours).
- 3. Prepare Citrate buffer fresh in the morning by combining 20 ml of 0.1 M Sodium citrate with 30 ml of 0.1 M Citrate acid to produce 0.1 M Citrate buffer. Adjust the pH to 4.0 by using 1 N NaOH. Filter-sterilize the Citrate buffer and store it in a sterile conical tube on ice.
- 4. Weigh the mice and ear tag them.
- 5. Add Citrate buffer to Streptozocin (STZ) in a sterile conical tube in a biosafety hood right before use to a final concentration of 20 mg/ml and keep on ice. For example, add 2.5 ml of Citrate buffer to 50 mg of Streptozocin.
- 6. Inject STZ at 200 mg/kg, i.p. (0.1 ml/10 g).
- 7. Inject Citrate buffer at 0.1 ml/10 g for the control groups.
- 8. Feed the mice with 10% sucrose in place of water on the 1st day of injection.
- 9. Remove the sucrose and feed the mice with fresh water the next day.
- 10. Monitor body weight, blood glucose and ketone bodies after injection as needed.
- 11. Harvest tissue samples after euthanize the mice with CO2.

Reagent and solution preparation:

Streptozocin powder

Weigh Streptozocin (Sigma 85882) and store at -20°C in aliquots (25-100 mg) in conical tubes that are sealed well with cap and parafilm, and clearly labeled with the chemical name and weight.

0.1 M sodium citrate

Dissolve 1.47 g of sodium citrate tribasic dihydrate (Sigma C0909, MW 294.10) in 50 ml ddH2O and store at room temperature.

0.1 M citric acid

Dissolve 1.05 g of citric acid monohydrate (Sigma C1909, MW 210.14) in 50 ml ddH2O and store at room temperature.