

VO2 treadmill protocol V1.3

This is a protocol for measurement of VO₂max in mice during treadmill running, which was modified from Columbus Instruments' manual and Lee-Young et al. (Wasserman) JBC, 2009; 284:23925-34.

Procedures

1. **Acclimate the mice** to the treadmill by running for 10 min at 10 m/min for 3 consecutive days at 9 am, prior to performing an analysis.
2. **Check Drierite** (should be blue). If purple, change it and SAVE purple Drierite for reuse.
3. **Turn on the power** (surge protector) closest to the tank to turn on the entire system.
4. **Turn on the computer** and login. Password is: yanlab.
5. **Weigh mice**. This factors into the VO₂ data since units are ml/kg/hr.
6. **Open CLAMS software**.
7. **Go to Setup Experiment** tab.
 - a. **Calibrate oxymax**
 - i. **Turn gas tank valve** to open. Pressure should be ~10-15 psi. DO NOT adjust the adjustment valve unless the pressure is not at the target pressure.
 - ii. **Input certified gas values** on the label of the gas tank. The values should be around 20.5% O₂ and 0.5% CO₂.
 - iii. **Click "Start Calibration"** in the CLAMS software when all lights of the equipment are stable (ready). This process will take approximately 8 minutes. You can set Treadmill Controller parameters during this time.
 - b. **Select the protocol** in "Treadmill Controller" - **See Below**.
 - c. **Select 3 cages** in the **Experiment setup** tab on the right side on the screen for the oxymax to sample gases. **Important**, if 4 cages are selected, it will rotate through all 4 cages.
 - d. **Click "Edit"** in the setup Experiment Tab.
 - i. **Click "Add subject"**. Enter a "Name" (Mouse number), "Mass" (weight), Strain (WT, AMPK KI, etc), Gender, Date of Birth
 - ii. **Click and drag** to a specific cage or treadmill channel
8. **Place mice** into designated treadmill channels. Make sure the latches (top and bottom) are secure and the treadmill channels are completely sealed.
9. **Turn on stimulators** at 0.3 mA and 1 Hz (lights will flash).
10. **Click "Data Log Location"**. Select location of file and name file.
11. **Click "start"** to start the experiment from **Setup Experiment or Run Experiment** tabs.
12. **Monitor the running activities** closely. **Turn off the electrical stimulator** when a mouse stays on the electrical grid for more than 5 seconds,.
13. **Click "stop"** 3 min after the last mouse stops running (allow for recording for all mice).
14. **Click "Export CSV"**.
 - a. Select cages used during the test
 - b. Select the tests/experiments to export
 - c. Save data file and excel files to: DESKTOP → CLAMS → DATE

VO_{2max} Protocol 3 Channel

Time (min)	Speed (m/min)	Incline (%)
-14	0	5
-12	5	5
0	10	5
3	13	5
7	16	5
10	19	5
13	22	5
16	25	5
19	28	5
22	31	5
25	34	5
28	37	5
31	40	5
34	43	5
37	46	5
40	49	5
43	52	5
46	55	5
49	58	5
52	61	5