Mouse voluntary running.v1

IMPORTANT!!

Never turn off the "Dataquest" program or the computer. You will disrupt other people's experiments!!! Please talk to me if you have any questions.

- 1. House mice in the Voluntary Running Room (Independence Park Facility Room 293) in a dark-light cycle (7:00 pm off, 7:00 am on) for 3 days in individual running cages with the running wheels locked by a straightened large paper clip.
- 2. Connect the cables between the matrix and the cages.
- 3. Log the information on the running wheel recording book.
- 4. Unlock the wheel around 5 pm on the 3rd day and set up the computer monitoring for each of the cages.
- 5. Double click on the "dog" signs of the designated channels, make sure the recording interval is 300 second (5') and then click "Ok".
- 6. Right click on the "dog" signs of the designated channels and click on "continuous sampling" to start the recording for the cages.
- 7. Unlock the wheel by removing the paper clip. Turn the wheel and wait for 5 min to see if the summary window shows a recorded number that is reasonable. For example, if you turn 5 times, the recording number should be 1 since 5 times/5 min = 5 turns/min.
- 8. Lock the wheels for 48 hours before sacrificing the mice if we want to determine the effect of long-term exercise training rather than acute effect of exercise.
- 9. Stop recording for each of the channels by right click on the "dog" sign and select "discontinue sampling". The information on the summary window for that channel will disappear.
- 10. After finish your study, move your data (only the data specific for your channels) to a new "Running Data" folder in a folder of your own. This will insure that other members could start to use the channels for their studies.

Calculations for Running Wheel Studies

- 1. Open Excel application
- 2. Open your file from Excel
- 3. Press "Fixed width" in dialog box
- 4. Press "Next" in dialog box (you should see columns)
- 5. Press "Finish" on dialog box to get to work sheet
- 6. Your data should be in columns (if not ask for help)
- 7. Multiply your raw data by 5 (this gives you rev/ 5 min)
- 8. Multiply rev/5 min by 0.357 (wheel circumference) (this gives you meters/5min)
- 9. Sum the meters/5min over 24 hour increments (this gives you meters/day)