Tamoxifen Injection

- 1. Add tamoxifen (Cat# T5648-1g, Sigma) at 10 mg/ml in peanut oil (P2144-250 ml, Sigma).
- 2. Incubate at 55°C water bath.
- 3. Votex every 5 min until completely dissolve.
- 4. Aliquot 0.5-1 ml into 1.5 ml tube.
- 5. Store at 4°C, which can be use for 3-4 days.
- 6. Use 23G needle to inject 100 μl (1 mg)/25 g of body weight (i.p.) once a day for 10 days.